

## BERNOULLI'S THEOREM APPARATUS

MAKE – LABTEK, ORIGIN – INDIA



## BERNOULLI'S THEOREM APPARATUS

### Supply Tank:

M.S. Powder Coating/M.S.FRP Lining/S.S. Unit of suitable size provided with piezo meter for the measurement of total potential head.

### Variable Cross Section Transparent Duct:

A transparent duct made of acrylic and of smooth variable cross section with suitably spaced 0.5m height replaceable high quality graduated glass tubes to measure the pressure head along the duct and a GM flow control valve.