

# Digital Muffle Furnance

MAKE – LABTEK, ORIGIN – INDIA



Light weight Outer body P.C.R.C. Sheet Duly Painted/Powder Coated, having Grooved refractory chamber to support Kanthal heating element on  $\frac{3}{4}$  sides as per requirement Temp. range  $1200^{\circ}\text{C}$ , working  $1150^{\circ}\text{C}$  fitted with high density ceramic fibre insulation to avoid heat/temp. loss, Fitted with Digital Temp Controller.

## OPTIONAL :

1. PID Temp. Controller can be fitted at an extra cost.