

## LEVEL CONTROL TRAINER MAKE – LABTEK, ORIGIN – INDIA



The Level Control System Trainer is the system which outlines the basics of closed loop level control process and other aspects related to it.

Some Additional Features: Compact Ergonomic Design. Computer Interface & SCADA software connectivity for analysis of Level Control System Trainer (Optional) System Dimensions: 4.5 Ft. (L) X 1.5 Ft. (W) X 4.5 Ft. (H) Services required: Water supply and drainage arrangement, Electric supply 230 V AC, 50 Hz. Clean, Dry Compressed Air Supply at 2.1 Kg/cm<sup>2</sup>.