

Water Bath-6 Holes

MAKE – LABTEK, ORIGIN – INDIA



Rectangular water bath double walled construction, inner reservoir of stainless steel. Outer body in mild steel with glass wool insulation in between, supplied with lid having one hole of 12.5cms dia & 6 holes of 7.5 cms dia. Concentric rings with lids are provided on holes. The top plate & concentric rings are of stainless steel temperature range ambient to 100° C. Controlled by capillary thermostat. Size of the chamber 30cm x 25cm x 12cm. These water baths are useful for laboratory work. Thermometers at extra cost suitable for 230 V A.C. Single Phase.

Customised sizes also available.