

CLOSED CIRCUIT FRANCIS TURBINE TEST RIG 5HP MAKE – LABTEK, ORIGIN – INDIA



Francis Turbine

Type	: Inward Flow Reaction Turbine
Capacity	: 5HP
Rated Speed	: 1250 RPM
Discharge Capacity	: 1800 Liters/Minute
Supply Head	: 25 Meters
Guide Vanes	: Gun Metal Vanes (Aerofoil Blade Shaped)
Loading	: Rope Brake
Material	: Cast Iron
Drum Size	: 300 mm diameter

Supply Pump Set