

CONSOLIDATION APP. (THREE GANG) MODEL NO – SL-GT-045, MAKE – LABTEK, ORIGIN – INDIA



IS 2720 (Part XV), IS 12287, BS 1377, ASTM D 2435

Rigidly manufactured from aluminous alloy casting to provide a high degree of accuracy with any frame distortions under load. The load bridge group is supported in high accuracy self-aligning seat balls. The beam provides three loading ratio: 9:1 10:1 11:1 and the beam assembly is fitted with an adjustable counterbalance weight. It consists of 3 consolidation test assemblies mounted on a single frame. Complete with 3 sets of Consolidation cell assembly for 60mm dia samples and three sets of weights each giving 10 kg/sq.cm. On each sample. And complete with rod holding the weights and coupling block holding the dial gauge or transducer. Supplied 'without': Consolidation cell, Weights, Dial Gauge (or transducer), holding bench which have to be ordered separately.

B-3A, Shiv Shakti Complex, East Vinod Nagar, 1st Floor, New Delhi - 110 091

Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915

E-mail: info@labtekindia.com, Web: www.labtekindia.com