

DIGITAL PHOTO COLORIMETER

MODEL SL-EE-001,002, MAKE – LABTEK, ORIGIN – INDIA



DIGITAL PHOTO COLORIMETER WITH AUTO ZERO CALIBRATION
(MODEL SL-EE-001)

DIGITAL PHOTO COLORIMETER WITH AUTO ZERO CALIBRATION (MODEL SL-EE-001)

STANDARD GLASS FILTER	8(Wavelength range 400-700nm).
MINIMUM VOLUME	1ml
DISPLAY	3 Digit Bright Red 7 Segment LED Display
KEY BOARD	1 key, Soft touch membrane type.
ABSORBANCE	0 to 1.99.
RESOLUTION	0.01.
ACCURACY	0.5% FSI.
DETECTOR	Selenium Photo Cell/Photodiode.
LIGHT SOURCE	6.3V-0.3Amp. Tungsten Lamp.
POWER	230 \hat{A} ± 10 % 50 Hz AC.

MICROPROCESSOR PHOTO COLORIMETER (MODEL SL-EE-002)