

DUAL OUTPUT DC POWER SUPPLY

MODEL NO – SL-EE-401, MAKE – LABTEK, ORIGIN – INDIA



Features: Should have the following features:

Dual independent adjustable DC outputs, Constant voltage and constant current, Independent, Series and parallel variety output, Low ripple and noise, Four LED displays with 3 digits, selectable for two sets output voltage and current display simultaneously, Auto tracking output

SPECIFICATIONS:

- **Line Regulation :** CV 1X10⁻⁴+3mV, CC 2x10⁻³+3mA
- **Load Regulation :** CV 1X10⁻⁴+4mV, CC 2x10⁻³+5mA
- **Ripple & Noise :** CV 1mVrms, CC 3mA_{rms}