

CLOSED CIRCUIT PELTON WHEEL TURBINE TEST RIG 5 HP MAKE – LABTEK, ORIGIN – INDIA



A. Pelton Wheel

Type	: Impulse
Capacity	: 5HP
Rated Speed	: 1000 RPM
Discharge Capacity	: 750 Liters/Minute
Supply Head	: 50 Meters
Buckets	: Gun Metal
B. Loading	: Rope Brake
Material	: Cast Iron
Drum Size	: 300 mm diameter

C. Supply Pump Set