

FLYWHEEL APPARATUS

MODEL NO – SL-MEL-35, MAKE – LABTEK, ORIGIN – INDIA



TECHNICAL SPECIFICATIONS :

Basic Wind Tunnel.

1. Self contained.
2. Wall mounting includes a pointer to the periphery of the wheel
3. To demonstrate the Second law of motion and energy storage
4. Running in ball bearings.
5. A load hanger and cord supplied.
6. 250mm x 30mm steel flywheel with engraved line on periphery
7. Instruction manual for student and lecturer provided.
8. Stop watch supplied
9. Set of weights.
10. Dimensions and Weight (excluding set of weights): 215 x 280 x 250 mm ; 19 kg (approx)