

HEAT TRANSFER LAB EQUIPMENTS THROUGH COMPOSITE WALL MAKE – LABTEK, ORIGIN – INDIA



Mild Steel, Asbestos & Brass SLabs

a. Thickness	:150 od 6 mm Thick
--------------	--------------------

02. Measuring Instruments

a. Digital Voltmeter Range	:0 to 300 V
----------------------------	-------------

b. Digital Ammeter & Range	:0 to 10 Amps
----------------------------	---------------

c. Digital Temperature Indicator with Selector switch