

## HYDRAULIC RAM TEST RIG CLOSED CIRCUIT HYDRAULIC RAM TEST RIG MAKE – LABTEK, ORIGIN – INDIA



### **Hydraulic Ram :**

Hydraulic Ram of size 50 x 15m to work against a delivery head up to 20m from a supply head of 2.5metres to discharge about 250LPH at 10m head.

### **Supply Tank :**

M.S Powder Coating/M.S FRP Lining/S.S Unit tank of suitable size with over flow arrangement for supply of water at constant head of 2.5 meters.



**Sun LabTek Equipments (I) Pvt. Ltd.**  
ISO:9001:2008 Certified Company

