

LOSS ON HEATING/ THIN FILM OVEN MODEL NO SL -BT -028, MAKE - LABTEK, ORIGIN - INDIA



Precise, Hot Air Drying Better Mineral / Blanket Insulation for high Temperature & to avoid heat loss Silent Hot Air Blower, Unique design of Air Circulation provide through out uniform Air movement Unique design of ventilation keeps the surface of the instruments from being Burnt even when the instruments i.e. Oven Temperature maintained at 200 ^\circ C Polish / Hair Line 304 grade S. Steel sheet interior, long operation, corrosion resistant Kanthal A-1 Super quality coil shaped & Air heater tubular model wound on side/on the back of the Oven for better accuracy Full feature with Digital Temperature Controller cum Indicator having Alarm facility (On Customer request) Toughened Glass view window to observe /Test the material without disturbing the Temperature condition of the chamber. Working Temperature required as per IS : ASTM is 163 $\text{^\circ C} \pm 1\text{^\circ C}$ Provide with detachable metal shaft (Both for Loss on Heating / Thin Film Oven) Reduction gear is fitted from outside rotated by a vertical shaft having 5-6 RPM.

B-3A, Shiv Shakti Complex, East Vinod Nagar, 1st Floor, New Delhi - 110 091

Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915

E-mail: info@labtekindia.com, Web: www.labtekindia.com