

MPFI-MULTICYLINDER-FOUR-STROKE-PETROL-ENGINE- TEST-RIG

MAKE – LABTEK, ORIGIN – INDIA



Introduction

The test rig is designed to provide self-contained facility for teaching Internal Combustion (spark Ignition) engine principles. The equipment is instrumented so that the following experiments could be performed.

- 1). Bhp Measurement
- 2). Ihp Measurement (By Morse Test Arrangement)
- 3). Fuel Consumption Measurement
- 4). Air Intake Measurement
- 5). Measurement Of Heat Rejected To Water Jacket
- 6). Heat Balance Test

B-3A, Shiv Shakti Complex, East Vinod Nagar, 1st Floor, New Delhi - 110 091

Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915

E-mail: info@labtekindia.com, Web: www.labtekindia.com