

ROTA METER TEST APPARATUS MAKE – LABTEK, ORIGIN – INDIA



Rota Meter:

Two Rota meters of size 20mm and 25mm for flow measurement fitted with respective pipeline

Measuring Tank:

M.S. Powder Coating/M.S.FRP Lining/S.S. Unit Measuring Tank of suitable size provided with overflow arrangement and a drain valve.

Sump Tank:

M.S. Powder Coating /MS.FRP Lining/SS. Unit of suitable size for suitable to store sufficient

B-3A, Shiv Shakti Complex, East Vinod Nagar, 1st Floor, New Delhi - 110 091

Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915

E-mail: info@labtekindia.com, Web: www.labtekindia.com