

SEMI AUTOMATIC COMPACTOR

MODEL NO – SL-GT-017, MAKE – LABTEK, ORIGIN – INDIA



IS 2720 (PART VII and VIII)

Now a days almost every soil laboratory conducts compaction and penetration tests on soils. In both standard and heavy compaction tests the operator has to count the number of blows of the rammer on each layer of soil during compaction. A mechanically cum electrically operated semi-automatic compactor is therefore devised to eliminate hand compaction process. Which also considerably saves time with the Rammers provided two stroke lengths, available, the compactor can be used for normal compaction/heavy compaction as well as C.B.R compaction.