


TORSION TESTING MACHINE MAKE – LABTEK, ORIGIN – INDIA



- 1). Unique "TALCON-CLAW" Grips for Round and Flat specimens
- 2). Torisonal Speed of 0.1/0.2 RPM
- 3). Torque Compensated Motors
- 4). Maintenance Free Gear Box
- 5). Auto Range Selection with Range Selector Wheel (Optional)
- 6). 3 or 4 Torque Ranges for Higher Accuracy (Optional)
- 7). Load-unload Capability
- 8). Variable Speed Drive (Optional)