


# CLOSED CIRCUIT APPARATUS FOR DETERMINATION OF FLOW THROUGH VENTURE METER AND ORIFICE METER

MAKE – LABTEK, ORIGIN – INDIA



### **Piping System:**

Consisting of Two G.I. Pipe lines, each line fitted with a venture meter & orifice meter of size 25 mm with a gate valve & provided with a set of distribution chamber having isolating cocks to facilitate experiment on individual venture meter or orifice meter.

### **Differential Head Measurement:**

Differential mercury manometer with mercury.

### **Measuring Tank:**

M.S. Powder Coating/M.S.FRP Lining/S.S. Unit of suitable size with an overflow arrangement & Provided with a PVC control valve & an elbow.